

311000 Mefix

Self-adhesive fabric

Product details

Size :	10 cm x 10m
Descriptive feature :	Self-adhesive fixation
Sterile :	Non-sterile
Brand :	Mefix

Images



Delivered items

311000-30

Sales released in: Algeria, Austria, Bahrain, Belgium, Brazil, Bulgaria, Canada, Chile, China, Croatia, Cyprus, Czechia, Denmark, Egypt, Estonia, Faroe Islands, Finland, France, Georgia, Germany, Gibraltar, Greece, Guadeloupe, Hong Kong, Hungary, Iceland, India, Iran (Islamic Republic of), Iraq, Ireland, Italy, Japan, Jordan, Kazakhstan, Kenya, Korea (the Republic of), Kuwait, Latvia, Lebanon, Lithuania, Luxembourg, Malaysia, Malta, Martinique, Morocco, Netherlands, New Zealand, Norway, Oman, Pakistan, Poland, Portugal, Qatar, Réunion, Romania, Saudi Arabia, Serbia, Singapore, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Taiwan (Province of China), Thailand, Turkey, Ukraine, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland, Viet Nam

Country of origin: Thailand

Shelf life: 3 years

Sterilization method: Non-sterile

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2023 Mölnlycke Health Care AB. All rights reserved.



transport box.

Is suitable for Tray: Yes

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430082374			
Consumer pack	1	7323190203352			
Dispenser box	1	7323190195459	74x106x70 mm		
Transport box	20	7323190195442	236x395x157 mm	14.6 dm3	3.6 / 3 kg
Pallet	2200	7323190195435	800x1200x1871 mm	1796.2 dm3	424.3 / 330.0 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Product Composition Film Products and Fixation Products

Product Component	Composition
Backing material	Hydrophobic polyester nonwoven
Adhesive	Polyacrylate
Protective release liner	Siliconised paper

Product Performance Fixation Products

Characteristics	Test Method	Internal Test Method	Unit	Requirement	Product Performance
Adhesion	ASTM D3330/D3330M-04, method F	T-1082	N/25 mm	Not specified	2.97
Moisture Vapour Transmission Rate (MVTR) of a wound dressing when in contact with water vapour	EN 13726-2 part 3:2	T-1070	g/m ² /24 h	Not specified	3448

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2023 Mölnlycke Health Care AB. All rights reserved.



Technical

Dimension

Dimension text	Dimension value
Product	10 cm x 10 m
Product	4 in x 11 yd

Classifications

Regulation type(s)	MDR Class I ns	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :	1			
Conformity Annexes :	IV			
Measuring Function :	No			
Notified body medical devices/PPE :	MPA	Local Authority Thailand	Local Authority EU Non EEC	Local Authority Thailand
510(k) clearance number :		NA		
Regulatory Released :	Austria, Belgium, Bulgaria, Croatia, Cyprus, Czechia, Denmark, Estonia, Faroe Islands, Finland, France, Germany, Gibraltar, Greece, Guadeloupe, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Martinique, Netherlands, Norway, Poland, Portugal, Réunion, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Thailand	Serbia	Thailand

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority India	Local Authority MEA	Local Authority Egypt	Local Authority MEA
510(k) clearance number :				NA
Regulatory Released :	India	Algeria, Bahrain, Iraq, Kenya, Lebanon, Morocco, Oman, South Africa, Turkey, Ukraine	Egypt	Turkey

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority South Korea	Local Authority Taiwan	Local Authority United Arab Emirates	Local Authority India
510(k) clearance number :				
Regulatory Released :	Korea (the Republic of)	Taiwan (Province of China)	United Arab Emirates	India

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority China	Local Authority MEA	Local Authority Pakistan	Local Authority Jordan
510(k) clearance number :		NA		
Regulatory Released :	China	Qatar	Pakistan	Jordan

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority Americas	Saudi Authority	Local Authority Kuwait	Local Authority India
510(k) clearance number :				
Regulatory Released :	Chile	Saudi Arabia	Kuwait	India

Find out more at www.molnlycke.com



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Class I
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority APAC	Local Authority MEA	Local Authority CIS	Health Canada
510(k) clearance number :		NA	NA	
Regulatory Released :	Hong Kong, Malaysia, New Zealand, Singapore	Algeria, Bahrain, Iraq, Kenya, Lebanon, Morocco, Oman, South Africa, Turkey, Ukraine	Kazakhstan	Canada

Classifications

Regulation type(s)	Class I	Non Medical Device
Intended Purpose :	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas	Mefix is a self-adhesive fabric for fixation of dressings, tubes and cannulas
MDR Classification Rule :		
Conformity Annexes :		
Measuring Function :		
Notified body medical devices/PPE :	Anvisa (Brazil)	PMDA (Japan)
510(k) clearance number :		
Regulatory Released :	Brazil	Japan

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN ISO 20417, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, EN ISO 10993-18, ISO 14001

Removable label

No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2023 Mölnlycke Health Care AB. All rights reserved.



GMDN Code (Global Medical Device Nomenclature)

16866 Basic adhesive tape

UNSPSC

42311552 Medical and surgical adherent tapes for general use

Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

CE Responsibility

Mölnlycke Health Care AB, Gamlestadvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000012JF

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2023 Mölnlycke Health Care AB. All rights reserved.

